



CRYSTAL

SPECIAL PASTE CERAMIC FOR MICRO-LAYERING ON
 LAYERED AND MONOLITHIC ZIRCONIA, LITHIUM SILICATE AND ALLOY-CERAMIC

- www.naturalceramic.eu - www.naturalceramic.it - www.naturalceramic.com - www.naturalceramic.it -



DENTAL CERAMIC



ITALIAN STYLE

Natural CRYSTAL

SPECIAL PASTE CERAMIC FOR MICRO-LAYERING

Natural CRYSTAL is the ceramic system in paste for micro-layering on any structure (alloy, zirconia, lithium silicate) compatible with the **Natural Ceramic System**.

Natural CRYSTAL is made of two customizable kits **CRYSTAL LIVE!** for the white aesthetic and **CRYSTAL GUM** for the pink aesthetic.

Natural CRYSTAL LIVE! is the ceramic paste system for the micro-layering of the cervical, dentin and incisal areas, the so-called "white area aesthetic".

Natural CRYSTAL GUM are ceramic masses in paste to recreate the aesthetics of soft tissues with the micro-stratification technique of the gingival areas, the so-called "pink area aesthetic".

Natural CRYSTAL LIVE!

FLUORESCENT CERAMIC PASTE FOR WHITE AESTHETIC

Natural CRYSTAL LIVE! is formed from:

- 4 neutral structure masses **LIVE! 1, LIVE! 2, LIVE! 3, LIVE! IWS** with different degrees of opacity;
- 2 special highly opaque masses for mamelons **LIVE! MAMELON CREAM** and **VANILLA**;
- 18 **LIVE! DENTIN** high fluorescence paste ceramic shades in A-D color range, each color from **A0** to **D4** available;
- 4 intensive masses for the incisal area **LIVE! ENAMEL S1, S2, S3, S4**;
- 6 transparent masses **LIVE! TRANSPA OPAL, CLEAR, PINK, BLUE, ORANGE** and **SICILY**.

LIVE! IWS

4 SPECIAL NEUTRAL CUSTOMIZABLE PASTE

LIVE! IWS ceramic paste masses can be used alone or mixed with Natural Stains and/or other CRYSTAL masses to desaturate them or vary their opacity



LIVE! MAMELON

2 SPECIAL MAMELON PASTE

LIVE! MAMELON are highly opaque paste shades to create effects within the micro-layering materials; these masses can also be used pure or mixed with the other Natural CRYSTAL paste and Natural



LIVE! DENTIN

18 SPECIAL A-D DENTIN PASTE

LIVE! DENTIN ceramic masses can be used to change the chromatic component of a zirconia restoration: characterizing, intensifying and/or modifying the basic shade of the element; LIVE! DENTIN pastes can also be mixed with Natural STAINS stains to obtain particularly sought-after effects.



LIVE! ENAMEL

4 SPECIAL ENAMEL PASTE

LIVE! ENAMEL are high visibility paste enamels, they emphasize the incisal area creating effects of depth even with minimal thicknesses.



LIVE! TRANSPARENT

6 SPECIAL TRANSPARENT paste

LIVE! TRANSPARENT are highly translucent and transparent masses, softer and more delicate, colored and opalescent, to create customizations in the incisal and also cervical areas.



Natural CRYSTAL GUM

SPECIAL PASTE CERAMIC FOR PINK AESTHETIC

Natural CRYSTAL GUM is formed from:

- 1 acromatic gingival mass **NEUTRAL** very useful to desaturate other gingival masses and/or make some white areas;
- 5 pink gingival masses **CHERRY, ORCHID, PEARL, SALMON** and **VELVET** with different degrees of saturation;
- 4 red gingival masses **ALFA, BURGUNDY, RUBY** and **WINE** with different degrees of color saturation, nuances and reflexes.

GUM NEUTRAL

1 SPECIAL GINGIVAL ACROMATIC PASTE

GUM NEUTRAL is an acromatic gingival mass for microlayering technique; bright, clear and with high value it is the best choice to light up other gingival masses.



GUM PINK

5 SPECIAL GINGIVAL PINK PASTE

GUM PINK are pink gingival masses, with several degrees of color saturation and reflections; they are ideal for replicating all the clearest tones of gum tissue.



GUM RED

4 SPECIAL GINGIVAL RED PASTE

GUM RED are gingival masses with a red chromatic base, with several degrees of saturation of the color and reflections; those ceramic pastes are ideal for reproducing all darker and vascularized tones of gingival tissue.



- www.naturalceramic.eu - www.naturalceramic.it - www.naturalceramic.com - www.naturalceramic.it -

 **TRESSIS ITALIA srl**

viale Italia, 194
31015 Conegliano (TV) ITALIA

Tel (+39) 0438 41 83 16

web www.tressis.it

email tressisitalia@gmail.com



DISTRIBUTED BY: